



Plasmon

Plasmon Library

New Diagnostic LEDs for
G Entry / Midrange

And

Enterprise G

RoHS Compliant Libraries



RoHS Libraries


Plasmon has added LEDs on certain circuit boards to enhance the serviceability of the Enterprise G and Gx libraries. The boards will only be available on RoHS compliant products. Plasmon libraries that are RoHS compliant will have a serial number of 20000 or greater.

The pages that follow describe the circuit boards with the new diagnostic LEDs.



G Entry / G Midrange

MTE Interface PCA PN: 910-100202-00




Board Silk Screen	LED Description	Blocked State
PA	Picker A: Home Position	On (home)
PB	Picker B: Home Position	On (home)
RI	Media Store Opto: Detects media in storage slot	On (media present)
AO	Auto Offset: Checks frame cutouts	On
MA	Media A: UDO VS MO Media Picker A	On (UDO media)
MB	Media B:UDO VS MO Media Picker B	On (UDO media)



G Entry / G Midrange

Lift PCA PN: 910-100203-00



Board Silk Screen	LED Description	Blocked State
FA	Flip A sensor blocked	On
FB	Flip B sensor blocked	On
LA	Lift Channel A Encoder	On/Off
LB	Lift Channel B Encoder	On/Off
WA	Flip / Pick Motor Channel A Encoder	On/Off
WB	Flip / Pick Motor Channel B Encoder	On/Off